

Reaction Tubes



Reaction Tubes

Overview: Max. Centrifuge Capacity in Technical Appendix


Further Coloured Caps p. 5 | 14

High chemical and temperature resistance

Available in different sizes

Available for Eppendorf, Vitatron, Roche and Beckman systems

DNase-free
RNase-free
human DNA-free
non-Pyrogenic



Cat.-No.	616 201	618 201	623 201	667 201	617 301	742 270
Description	reaction tube	reaction tube	reaction tube	reaction tube	reaction tube	"Cobas" cup
Colour	natural	natural	natural	natural	natural	blue
Volume [ml]	1.5	1.5	2.0	0.5	0.4	0.7
Graduation	+	-	+	-	-	-
Cap, attached	+	-	+	+	+	+
Suitable for	Eppendorf	Eppendorf	Eppendorf	Vitatron	Beckman	Roche
Sterile	-	-	-	-	-	-
Quantity per bag/case	500/4000	500/3000	500/4000	1000/10000	1000/5000	500/15000



Cat.-No.	716 201	717 201	722 201
Description	reaction tube	reaction tube	reaction tube
Volume [ml]	1.5	1.5	2.0
Graduation	-	-	-
Screw cap, Cat.-No.	366 380	367 XXX	366 380
Skirt	-	+	+
Sterile	-	-	-
Quantity per bag/case	500/5000	500/5000	500/5000

Michał Hadryś
+ 48 695 222 099
michal.hadrys@panalytica.pl